

Antibodies Product Datasheet

Product Name:	MAD2L2 Antibody
Catalog Number:	C48134
Reactivity:	Human, Mouse, Rat
Applications:	WB
Dilutions:	WB 1:500-2000
Modification:	Unmodified/Total
Source:	Polyclonal, Rabbit,IgG
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Immunogen:	Synthesized peptide derived from human MAD2L2
Concentration:	1 mg/ml
Gene Name:	MAD2L2 MAD2B REV7
Protein Name:	Mitotic spindle assembly checkpoint protein MAD2B (Mitotic arrest deficient 2-like protein 2) (MAD2-like protein 2) (REV7 homolog) (hREV7)
Other Name:	
Human Swiss Prot No.:	Q9UI95
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	This antibody detects endogenous levels of MAD2L2 at Human, Mouse,Rat
Storage and Stability:	Store at -20° C/1 year