

Antibodies Product Datasheet

Product Name:	CCP1 Antibody
Catalog Number:	C48135
Reactivity:	Human, Mouse
Applications:	WB
Dilutions:	WB 1:500-2000
Modification:	Unmodified/Total
Source:	Polyclonal, Rabbit,IgG
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Immunogen:	Synthesized peptide derived from human CCP1
Concentration:	1 mg/ml
Gene Name:	AGTPBP1 CCP1 KIAA1035 NNA1
Protein Name:	Cytosolic carboxypeptidase 1 (EC 3.4.17.-) (ATP/GTP-binding protein 1) (Nervous system nuclear protein induced by axotomy protein 1 homolog)
Other Name:	
Human Swiss Prot No.:	Q9UPW5
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	This antibody detects endogenous levels of CCP1 at Human, Mouse
Storage and Stability:	Store at -20° C/1 year