

Antibodies Product Datasheet

Product Name:	LMOD1 Antibody
Catalog Number:	C48136
Reactivity:	Human, Mouse
Applications:	WB
Dilutions:	WB 1:500-2000
Modification:	Unmodified/Total
Source:	Polyclonal, Rabbit, IgG
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Immunogen:	Synthesized peptide derived from human LMOD1
Concentration:	1 mg/ml
Gene Name:	LMOD1
Protein Name:	Leiomodin-1 (64 kDa autoantigen 1D) (64 kDa autoantigen 1D3) (64 kDa autoantigen D1) (Leiomodin, muscle form) (Smooth muscle leiomodin) (SM-Lmod) (Thyroid-associated ophthalmopathy autoantigen)
Other Name:	
Human Swiss Prot No.:	P29536
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	This antibody detects endogenous levels of LMOD1 at Human, Mouse
Storage and Stability:	Store at -20° C/1 year