

Antibodies Product Datasheet

Product Name:	L2HGDH Antibody
Catalog Number:	C48137
Reactivity:	Human, Mouse
Applications:	WB
Dilutions:	WB 1:500-2000
Modification:	Unmodified/Total
Source:	Polyclonal, Rabbit,IgG
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Immunogen:	Synthesized peptide derived from human L2HGDH
Concentration:	1 mg/ml
Gene Name:	L2HGDH C14orf160
Protein Name:	L-2-hydroxyglutarate dehydrogenase, mitochondrial (EC 1.1.99.2) (Duranin)
Other Name:	
Human Swiss Prot No.:	Q9H9P8
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	This antibody detects endogenous levels of L2HGDH at Human, Mouse
Storage and Stability:	Store at -20° C/1 year