

Antibodies Product Datasheet

Product Name:	Choline Kinase Alpha Antibody
Catalog Number:	C48142
Reactivity:	Human, Mouse, Rat
Applications:	WB
Dilutions:	WB 1:500-2000
Modification:	Unmodified/Total
Source:	Polyclonal, Rabbit,IgG
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Immunogen:	Synthesized peptide derived from human Choline Kinase Alpha
Concentration:	1 mg/ml
Gene Name:	CHKA CHK CKI
Protein Name:	Choline kinase alpha (CK) (EC 2.7.1.32) (CHETK-alpha) (Ethanolamine kinase) (EK) (EC 2.7.1.82)
Other Name:	
Human Swiss Prot No.:	P35790
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	This antibody detects endogenous levels of Choline Kinase Alpha at Human, Mouse,Rat
Storage and Stability:	Store at -20° C/1 year