

Antibodies Product Datasheet

Product Name:	GEFT Antibody
Catalog Number:	C48146
Reactivity:	Human, Mouse, Rat
Applications:	WB
Dilutions:	WB 1:500-2000
Modification:	Unmodified/Total
Source:	Polyclonal, Rabbit,IgG
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Immunogen:	Synthesized peptide derived from human GEFT
Concentration:	1 mg/ml
Gene Name:	ARHGEF25 GEFT
Protein Name:	Rho guanine nucleotide exchange factor 25 (Guanine nucleotide exchange factor GEFT) (Rac/Cdc42/Rho exchange factor GEFT) (RhoA/Rac/Cdc42 guanine nucleotide exchange factor GEFT) (p63RhoGEF)
Other Name:	
Human Swiss Prot No.:	Q86VW2
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	This antibody detects endogenous levels of GEFT at Human, Mouse,Rat
Storage and Stability:	Store at -20° C/1 year