

Antibodies Product Datasheet

Product Name:	BCAT2 Antibody
Catalog Number:	C48162
Reactivity:	Human, Mouse, Rat
Applications:	WB
Dilutions:	WB 1:500-2000
Modification:	Unmodified/Total
Source:	Polyclonal, Rabbit,IgG
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Immunogen:	Synthesized peptide derived from human BCAT2
Concentration:	1 mg/ml
Gene Name:	BCAT2 BCATM BCT2 ECA40
Protein Name:	Branched-chain-amino-acid aminotransferase, mitochondrial (BCAT(m)) (EC 2.6.1.42) (Placental protein 18) (PP18)
Other Name:	
Human Swiss Prot No.:	O15382
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	This antibody detects endogenous levels of BCAT2 at Human, Mouse,Rat
Storage and Stability:	Store at -20° C/1 year