

Antibodies Product Datasheet

Product Name:	ASGR1 Antibody
Catalog Number:	C48171
Reactivity:	Human, Mouse, Rat
Applications:	WB
Dilutions:	WB 1:500-2000
Modification:	Unmodified/Total
Source:	Polyclonal, Rabbit,IgG
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Immunogen:	Synthesized peptide derived from human ASGR1
Concentration:	1 mg/ml
Gene Name:	ASGR1 CLEC4H1
Protein Name:	Asialoglycoprotein receptor 1 (ASGP-R 1) (ASGPR 1) (C-type lectin domain family 4 member H1) (Hepatic lectin H1) (HL-1)
Other Name:	
Human Swiss Prot No.:	P07306
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	This antibody detects endogenous levels of ASGR1 at Human, Mouse,Rat
Storage and Stability:	Store at -20° C/1 year