

Antibodies Product Datasheet

Product Name:	ACAA1 Antibody
Catalog Number:	C48184
Reactivity:	Human
Applications:	WB
Dilutions:	WB 1:500-2000
Modification:	Unmodified/Total
Source:	Polyclonal, Rabbit,IgG
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Immunogen:	Synthesized peptide derived from human ACAA1
Concentration:	1 mg/ml
Gene Name:	ACAA1 ACAA PTHIO
Protein Name:	3-ketoacyl-CoA thiolase, peroxisomal (EC 2.3.1.16) (Acetyl-CoA acyltransferase) (Beta-ketothiolase) (Peroxisomal 3-oxoacyl-CoA thiolase)
Other Name:	
Human Swiss Prot No.:	P09110
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	This antibody detects endogenous levels of ACAA1 at Human
Storage and Stability:	Store at -20° C/1 year