

Antibodies Product Datasheet

Product Name:	HLA-F Antibody
Catalog Number:	C48191
Reactivity:	Human
Applications:	WB
Dilutions:	WB 1:500-2000
Modification:	Unmodified/Total
Source:	Polyclonal, Rabbit,IgG
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Immunogen:	Synthesized peptide derived from human HLA-F
Concentration:	1 mg/ml
Gene Name:	HLA-F HLA-5.4 HLA-F
Protein Name:	HLA class I histocompatibility antigen, alpha chain F (CDA12) (HLA F antigen) (Leukocyte antigen F) (MHC class I antigen F)
Other Name:	
Human Swiss Prot No.:	P30511
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	This antibody detects endogenous levels of HLA-F at Human
Storage and Stability:	Store at -20° C/1 year