

Antibodies Product Datasheet

Product Name:	KATNA1 Antibody
Catalog Number:	C48203
Reactivity:	Human, Rat
Applications:	WB
Dilutions:	WB 1:500-2000
Modification:	Unmodified/Total
Source:	Polyclonal, Rabbit,IgG
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Immunogen:	Synthesized peptide derived from human KATNA1
Concentration:	1 mg/ml
Gene Name:	KATNA1
Protein Name:	Katanin p60 ATPase-containing subunit A1 (Katanin p60 subunit A1) (EC 3.6.4.3) (p60 katanin)
Other Name:	
Human Swiss Prot No.:	O75449
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	This antibody detects endogenous levels of KATNA1 at Human,Rat
Storage and Stability:	Store at -20° C/1 year