

Antibodies Product Datasheet

Product Name:	CPT1C Antibody
Catalog Number:	C48204
Reactivity:	Human, Mouse
Applications:	WB
Dilutions:	WB 1:500-2000
Modification:	Unmodified/Total
Source:	Polyclonal, Rabbit,IgG
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Immunogen:	Synthesized peptide derived from human CPT1C
Concentration:	1 mg/ml
Gene Name:	CPT1C CATL1
Protein Name:	Carnitine O-palmitoyltransferase 1, brain isoform (CPT1-B) (EC 2.3.1.21) (CPT 1C) (Carnitine O-palmitoyltransferase I, brain isoform) (CPTI-B) (Carnitine palmitoyltransferase 1C)
Other Name:	
Human Swiss Prot No.:	Q8TCG5
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	This antibody detects endogenous levels of CPT1C at Human, Mouse
Storage and Stability:	Store at -20° C/1 year