

Antibodies Product Datasheet

Product Name:	HAX1 Antibody
Catalog Number:	C48207
Reactivity:	Human, Mouse, Rat
Applications:	WB
Dilutions:	WB 1:500-2000
Modification:	Unmodified/Total
Source:	Polyclonal, Rabbit,IgG
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Immunogen:	Synthesized peptide derived from human HAX1
Concentration:	1 mg/ml
Gene Name:	HAX1 HS1BP1
Protein Name:	HCLS1-associated protein X-1 (HS1-associating protein X-1) (HAX-1) (HS1-binding protein 1) (HSP1BP-1)
Other Name:	
Human Swiss Prot No.:	O00165
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	This antibody detects endogenous levels of HAX1 at Human, Mouse,Rat
Storage and Stability:	Store at -20° C/1 year