

Antibodies Product Datasheet

Product Name:	DHX36 Antibody
Catalog Number:	C48208
Reactivity:	Human, Mouse
Applications:	WB
Dilutions:	WB 1:500-2000
Modification:	Unmodified/Total
Source:	Polyclonal, Rabbit, IgG
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Immunogen:	Synthesized peptide derived from human DHX36
Concentration:	1 mg/ml
Gene Name:	DHX36 DDX36 KIAA1488 MLEL1 RHAU
Protein Name:	Probable ATP-dependent RNA helicase DHX36 (EC 3.6.4.13) (DEAH box protein 36) (MLE-like protein 1) (RNA helicase associated with AU-rich element ARE)
Other Name:	
Human Swiss Prot No.:	Q9H2U1
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	This antibody detects endogenous levels of DHX36 at Human, Mouse
Storage and Stability:	Store at -20° C/1 year