

Antibodies Product Datasheet

Product Name:	DSS1 Antibody
Catalog Number:	C48218
Reactivity:	Human, Mouse
Applications:	WB
Dilutions:	WB 1:500-2000
Modification:	Unmodified/Total
Source:	Polyclonal, Rabbit,IgG
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Immunogen:	Synthesized peptide derived from human DSS1
Concentration:	1 mg/ml
Gene Name:	SHFM1 DSS1 SHFDG1
Protein Name:	26S proteasome complex subunit DSS1 (Deleted in split hand/split foot protein 1) (Split hand/foot deleted protein 1) (Split hand/foot malformation type 1 protein)
Other Name:	
Human Swiss Prot No.:	P60896
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	This antibody detects endogenous levels of DSS1 at Human, Mouse
Storage and Stability:	Store at -20° C/1 year