

Antibodies Product Datasheet

Product Name:	ETFDH Antibody
Catalog Number:	C48227
Reactivity:	Human, Mouse, Rat
Applications:	WB
Dilutions:	WB 1:500-2000
Modification:	Unmodified/Total
Source:	Polyclonal, Rabbit,IgG
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Immunogen:	Synthesized peptide derived from human ETFDH
Concentration:	1 mg/ml
Gene Name:	ETFDH
Protein Name:	Electron transfer flavoprotein-ubiquinone oxidoreductase, mitochondrial (ETF-QO) (ETF-ubiquinone oxidoreductase) (EC 1.5.5.1) (Electron-transferring-flavoprotein dehydrogenase) (ETF dehydrogenase)
Other Name:	
Human Swiss Prot No.:	Q16134
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	This antibody detects endogenous levels of ETFDH at Human, Mouse,Rat
Storage and Stability:	Store at -20° C/1 year