

Antibodies Product Datasheet

Product Name:	DHCR24 Antibody
Catalog Number:	C48239
Reactivity:	Human, Mouse, Rat
Applications:	WB
Dilutions:	WB 1:500-2000
Modification:	Unmodified/Total
Source:	Polyclonal, Rabbit,IgG
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Immunogen:	Synthesized peptide derived from human DHCR24
Concentration:	1 mg/ml
Gene Name:	DHCR24 KIAA0018
Protein Name:	Delta(24)-sterol reductase (EC 1.3.1.72) (24-dehydrocholesterol reductase) (3-beta-hydroxysterol delta-24-reductase) (Diminuto/dwarf1 homolog) (Seladin-1)
Other Name:	
Human Swiss Prot No.:	Q15392
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	This antibody detects endogenous levels of DHCR24 at Human, Mouse,Rat
Storage and Stability:	Store at -20° C/1 year