

Antibodies Product Datasheet

Product Name:	IPO4 Antibody
Catalog Number:	C48248
Reactivity:	Human, Mouse
Applications:	WB
Dilutions:	WB 1:500-2000
Modification:	Unmodified/Total
Source:	Polyclonal, Rabbit,IgG
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Immunogen:	Synthesized peptide derived from human IPO4
Concentration:	1 mg/ml
Gene Name:	IPO4 IMP4B RANBP4
Protein Name:	Importin-4 (Imp4) (Importin-4b) (Imp4b) (Ran-binding protein 4) (RanBP4)
Other Name:	
Human Swiss Prot No.:	Q8TEX9
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	This antibody detects endogenous levels of IPO4 at Human, Mouse
Storage and Stability:	Store at -20° C/1 year