

Antibodies Product Datasheet

Product Name:	ENOX2 Antibody
Catalog Number:	C48251
Reactivity:	Human, Mouse
Applications:	WB
Dilutions:	WB 1:500-2000
Modification:	Unmodified/Total
Source:	Polyclonal, Rabbit,IgG
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Immunogen:	Synthesized peptide derived from human ENOX2
Concentration:	1 mg/ml
Gene Name:	ENOX2 COVA1
Protein Name:	Ecto-NOX disulfide-thiol exchanger 2 (APK1 antigen) (Cytosolic ovarian carcinoma antigen 1) (Tumor-associated hydroquinone oxidase) (tNOX) [Includes: Hydroquinone [NADH] oxidase (EC 1.-.-.)Protein disulfide-thiol oxidoreductase (EC 1.-.-.)]
Other Name:	
Human Swiss Prot No.:	Q16206
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	This antibody detects endogenous levels of ENOX2 at Human, Mouse
Storage and Stability:	Store at -20° C/1 year