

Antibodies Product Datasheet

Product Name:	HNRNPUL1 Antibody
Catalog Number:	C48266
Reactivity:	Human, Mouse
Applications:	WB
Dilutions:	WB 1:500-2000
Modification:	Unmodified/Total
Source:	Polyclonal, Rabbit,IgG
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Immunogen:	Synthesized peptide derived from human HNRNPUL1
Concentration:	1 mg/ml
Gene Name:	HNRNPUL1 E1BAP5 HNRPUL1
Protein Name:	Heterogeneous nuclear ribonucleoprotein U-like protein 1 (Adenovirus early region 1B-associated protein 5) (E1B-55 kDa-associated protein 5) (E1B-AP5)
Other Name:	
Human Swiss Prot No.:	Q9BUJ2
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	This antibody detects endogenous levels of HNRNPUL1 at Human, Mouse
Storage and Stability:	Store at -20° C/1 year