

Antibodies Product Datasheet

Product Name:	GP73/GOLPH2 Antibody
Catalog Number:	C48273
Reactivity:	Human, Mouse
Applications:	WB
Dilutions:	WB 1:500-2000
Modification:	Unmodified/Total
Source:	Polyclonal, Rabbit,IgG
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Immunogen:	Synthesized peptide derived from human GP73/GOLPH2
Concentration:	1 mg/ml
Gene Name:	GOLM1 C9orf155 GOLPH2 PSEC0242 UNQ686/PRO1326
Protein Name:	Golgi membrane protein 1 (Golgi membrane protein GP73) (Golgi phosphoprotein 2)
Other Name:	
Human Swiss Prot No.:	Q8NBJ4
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	This antibody detects endogenous levels of GP73/GOLPH2 at Human, Mouse
Storage and Stability:	Store at -20° C/1 year