

Antibodies Product Datasheet

Product Name:	KCTD12 Antibody
Catalog Number:	C48302
Reactivity:	Human, Mouse
Applications:	WB
Dilutions:	WB 1:500-2000
Modification:	Unmodified/Total
Source:	Polyclonal, Rabbit,IgG
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Immunogen:	Synthesized peptide derived from human KCTD12
Concentration:	1 mg/ml
Gene Name:	KCTD12 C13orf2 KIAA1778 PFET1
Protein Name:	BTB/POZ domain-containing protein KCTD12 (Pfetin) (Predominantly fetal expressed T1 domain)
Other Name:	
Human Swiss Prot No.:	Q96CX2
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	This antibody detects endogenous levels of KCTD12 at Human, Mouse
Storage and Stability:	Store at -20° C/1 year