

Antibodies Product Datasheet

Product Name:	PSMB1 Antibody
Catalog Number:	C48305
Reactivity:	Human, Mouse, Rat
Applications:	WB
Dilutions:	WB 1:500-2000
Modification:	Unmodified/Total
Source:	Polyclonal, Rabbit,IgG
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Immunogen:	Synthesized peptide derived from human PSMB1
Concentration:	1 mg/ml
Gene Name:	PSMB1 PSC5
Protein Name:	Proteasome subunit beta type-1 (EC 3.4.25.1) (Macropain subunit C5) (Multicatalytic endopeptidase complex subunit C5) (Proteasome component C5) (Proteasome gamma chain)
Other Name:	
Human Swiss Prot No.:	P20618
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	This antibody detects endogenous levels of PSMB1 at Human, Mouse,Rat
Storage and Stability:	Store at -20° C/1 year