

Antibodies Product Datasheet

Product Name:	DDX39B Antibody
Catalog Number:	C48316
Reactivity:	Human, Mouse, Rat
Applications:	WB
Dilutions:	WB 1:500-2000
Modification:	Unmodified/Total
Source:	Polyclonal, Rabbit,IgG
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Immunogen:	Synthesized peptide derived from human DDX39B
Concentration:	1 mg/ml
Gene Name:	DDX39B BAT1 UAP56
Protein Name:	Spliceosome RNA helicase DDX39B (EC 3.6.4.13) (56 kDa U2AF65-associated protein) (ATP-dependent RNA helicase p47) (DEAD box protein UAP56) (HLA-B-associated transcript 1 protein)
Other Name:	
Human Swiss Prot No.:	Q13838
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	This antibody detects endogenous levels of DDX39B at Human, Mouse,Rat
Storage and Stability:	Store at -20° C/1 year