

Antibodies Product Datasheet

Product Name:	ADRM1 Antibody
Catalog Number:	C48319
Reactivity:	Human, Mouse, Rat
Applications:	WB
Dilutions:	WB 1:500-2000
Modification:	Unmodified/Total
Source:	Polyclonal, Rabbit,IgG
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Immunogen:	Synthesized peptide derived from human ADRM1
Concentration:	1 mg/ml
Gene Name:	ADRM1 GP110
Protein Name:	Proteasomal ubiquitin receptor ADRM1 (110 kDa cell membrane glycoprotein) (Gp110) (Adhesion-regulating molecule 1) (ARM-1) (Proteasome regulatory particle non-ATPase 13) (hRpn13) (Rpn13 homolog)
Other Name:	
Human Swiss Prot No.:	Q16186
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	This antibody detects endogenous levels of ADRM1 at Human, Mouse,Rat
Storage and Stability:	Store at -20° C/1 year