

Antibodies Product Datasheet

Product Name:	BLOC1S1 Antibody
Catalog Number:	C48322
Reactivity:	Human, Mouse, Rat
Applications:	WB
Dilutions:	WB 1:500-2000
Modification:	Unmodified/Total
Source:	Polyclonal, Rabbit,IgG
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Immunogen:	Synthesized peptide derived from human BLOC1S1
Concentration:	1 mg/ml
Gene Name:	BLOC1S1 BLOS1 GCN5L1 RT14
Protein Name:	Biogenesis of lysosome-related organelles complex 1 subunit 1 (BLOC-1 subunit 1) (GCN5-like protein 1) (Protein RT14)
Other Name:	
Human Swiss Prot No.:	P78537
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	This antibody detects endogenous levels of BLOC1S1 at Human, Mouse,Rat
Storage and Stability:	Store at -20° C/1 year