

Antibodies Product Datasheet

Product Name:	FBXW7 Antibody
Catalog Number:	C48331
Reactivity:	Human, Mouse
Applications:	WB
Dilutions:	WB 1:500-2000
Modification:	Unmodified/Total
Source:	Polyclonal, Rabbit, IgG
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Immunogen:	Synthesized peptide derived from human FBXW7
Concentration:	1 mg/ml
Gene Name:	FBXW7 FBW7 FBX30 SEL10
Protein Name:	F-box/WD repeat-containing protein 7 (Archipelago homolog) (hAgo) (F-box and WD-40 domain-containing protein 7) (F-box protein FBX30) (SEL-10) (hCdc4)
Other Name:	
Human Swiss Prot No.:	Q969H0
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	This antibody detects endogenous levels of FBXW7 at Human, Mouse
Storage and Stability:	Store at -20° C/1 year