

Antibodies Product Datasheet

Product Name:	GART Antibody
Catalog Number:	C48335
Reactivity:	Human, Mouse
Applications:	WB
Dilutions:	WB 1:500-2000
Modification:	Unmodified/Total
Source:	Polyclonal, Rabbit,IgG
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Immunogen:	Synthesized peptide derived from human GART
Concentration:	1 mg/ml
Gene Name:	GART PGFT PRGS
Protein Name:	Trifunctional purine biosynthetic protein adenosine-3 [Includes: Phosphoribosylamine--glycine ligase (EC 6.3.4.13) (Glycinamide ribonucleotide synthetase) (GARS) (Phosphoribosylglycinamide synthetase)Phosphoribosylformylglycinamidine cyclo-ligase (EC 6.3.
Other Name:	
Human Swiss Prot No.:	P22102
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	This antibody detects endogenous levels of GART at Human, Mouse
Storage and Stability:	Store at -20° C/1 year