

Antibodies Product Datasheet

Product Name:	GPR177 Antibody
Catalog Number:	C48337
Reactivity:	Human, Mouse, Rat
Applications:	WB
Dilutions:	WB 1:500-2000
Modification:	Unmodified/Total
Source:	Polyclonal, Rabbit,IgG
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Immunogen:	Synthesized peptide derived from human GPR177
Concentration:	1 mg/ml
Gene Name:	WLS C1orf139 GPR177 UNQ85/PRO18667
Protein Name:	Protein wntless homolog (Integral membrane protein GPR177) (Protein evenness interrupted homolog) (EVI) (Putative NF-kappa-B-activating protein 373)
Other Name:	
Human Swiss Prot No.:	Q5T9L3
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	This antibody detects endogenous levels of GPR177 at Human, Mouse,Rat
Storage and Stability:	Store at -20° C/1 year