

Antibodies Product Datasheet

Product Name:	ATP6V0A1 Antibody
Catalog Number:	C48338
Reactivity:	Human, Mouse, Rat
Applications:	WB
Dilutions:	WB 1:500-2000
Modification:	Unmodified/Total
Source:	Polyclonal, Rabbit,IgG
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Immunogen:	Synthesized peptide derived from human ATP6V0A1
Concentration:	1 mg/ml
Gene Name:	ATP6V0A1 ATP6N1 ATP6N1A VPP1
Protein Name:	V-type proton ATPase 116 kDa subunit a isoform 1 (V-ATPase 116 kDa isoform a1) (Clathrin-coated vesicle/synaptic vesicle proton pump 116 kDa subunit) (Vacuolar adenosine triphosphatase subunit Ac116) (Vacuolar proton pump subunit 1) (Vacuolar proton trans
Other Name:	
Human Swiss Prot No.:	Q93050
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	This antibody detects endogenous levels of ATP6V0A1 at Human, Mouse,Rat
Storage and Stability:	Store at -20° C/1 year