

Antibodies Product Datasheet

Product Name:	KLF6 Antibody
Catalog Number:	C48350
Reactivity:	Human, Rat
Applications:	WB
Dilutions:	WB 1:500-2000
Modification:	Unmodified/Total
Source:	Polyclonal, Rabbit,IgG
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Immunogen:	Synthesized peptide derived from human KLF6
Concentration:	1 mg/ml
Gene Name:	KLF6 BCD1 COPEB CPBP ST12
Protein Name:	Krueppel-like factor 6 (B-cell-derived protein 1) (Core promoter element-binding protein) (GC-rich sites-binding factor GBF) (Proto-oncogene BCD1) (Suppressor of tumorigenicity 12 protein) (Transcription factor Zf9)
Other Name:	
Human Swiss Prot No.:	Q99612
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	This antibody detects endogenous levels of KLF6 at Human,Rat
Storage and Stability:	Store at -20° C/1 year