

Antibodies Product Datasheet

Product Name:	GPR137 Antibody
Catalog Number:	C48351
Reactivity:	Human, Mouse
Applications:	WB
Dilutions:	WB 1:500-2000
Modification:	Unmodified/Total
Source:	Polyclonal, Rabbit,IgG
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Immunogen:	Synthesized peptide derived from human GPR137
Concentration:	1 mg/ml
Gene Name:	GPR137 C11orf4 GPR137A TM7SF1L1
Protein Name:	Integral membrane protein GPR137 (Transmembrane 7 superfamily member 1-like 1 protein)
Other Name:	
Human Swiss Prot No.:	Q96N19
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	This antibody detects endogenous levels of GPR137 at Human, Mouse
Storage and Stability:	Store at -20° C/1 year