

Antibodies Product Datasheet

Product Name:	ATAD3A/B Antibody
Catalog Number:	C48362
Reactivity:	Human
Applications:	WB
Dilutions:	WB 1:500-2000
Modification:	Unmodified/Total
Source:	Polyclonal, Rabbit, IgG
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Immunogen:	Synthesized peptide derived from human ATAD3A/B
Concentration:	1 mg/ml
Gene Name:	ATAD3B KIAA1273 TOB3
Protein Name:	ATPase family AAA domain-containing protein 3B (AAA-TOB3)
Other Name:	
Human Swiss Prot No.:	Q5T9A4
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	This antibody detects endogenous levels of ATAD3A/B at Human
Storage and Stability:	Store at -20° C/1 year