

Antibodies Product Datasheet

Product Name:	LAPTM4B Antibody
Catalog Number:	C48367
Reactivity:	Human, Mouse, Rat
Applications:	WB
Dilutions:	WB 1:500-2000
Modification:	Unmodified/Total
Source:	Polyclonal, Rabbit,IgG
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Immunogen:	Synthesized peptide derived from human LAPTM4B
Concentration:	1 mg/ml
Gene Name:	LAPTM4B PSEC0001
Protein Name:	Lysosomal-associated transmembrane protein 4B (Lysosome-associated transmembrane protein 4-beta)
Other Name:	
Human Swiss Prot No.:	Q86VI4
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	This antibody detects endogenous levels of LAPTM4B at Human, Mouse,Rat
Storage and Stability:	Store at -20° C/1 year