

Antibodies Product Datasheet

Product Name:	ARHGAP10 Antibody
Catalog Number:	C48373
Reactivity:	Human, Mouse
Applications:	WB
Dilutions:	WB 1:500-2000
Modification:	Unmodified/Total
Source:	Polyclonal, Rabbit,IgG
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Immunogen:	Synthesized peptide derived from human ARHGAP10
Concentration:	1 mg/ml
Gene Name:	ARHGAP10 GRAF2
Protein Name:	Rho GTPase-activating protein 10 (GTPase regulator associated with focal adhesion kinase 2) (Graf-related protein 2) (Rho-type GTPase-activating protein 10)
Other Name:	
Human Swiss Prot No.:	A1A4S6
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	This antibody detects endogenous levels of ARHGAP10 at Human, Mouse
Storage and Stability:	Store at -20° C/1 year