

Antibodies Product Datasheet

Product Name:	GPT/ ALT1 Antibody
Catalog Number:	C48378
Reactivity:	Human, Mouse, Rat
Applications:	WB
Dilutions:	WB 1:500-2000
Modification:	Unmodified/Total
Source:	Polyclonal, Rabbit,IgG
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Immunogen:	Synthesized peptide derived from human GPT/ ALT1
Concentration:	1 mg/ml
Gene Name:	GPT AAT1 GPT1
Protein Name:	Alanine aminotransferase 1 (ALT1) (EC 2.6.1.2) (Glutamate pyruvate transaminase 1) (GPT 1) (Glutamic--alanine transaminase 1) (Glutamic--pyruvic transaminase 1)
Other Name:	
Human Swiss Prot No.:	P24298
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	This antibody detects endogenous levels of GPT/ ALT1 at Human, Mouse,Rat
Storage and Stability:	Store at -20° C/1 year