

Antibodies Product Datasheet

Product Name:	GABPB1 Antibody
Catalog Number:	C48387
Reactivity:	Human, Mouse
Applications:	WB
Dilutions:	WB 1:500-2000
Modification:	Unmodified/Total
Source:	Polyclonal, Rabbit,IgG
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Immunogen:	Synthesized peptide derived from human GABPB1
Concentration:	1 mg/ml
Gene Name:	GABPB1 E4TF1B GABPB GABPB2
Protein Name:	GA-binding protein subunit beta-1 (GABP subunit beta-1) (GABPB-1) (GABP subunit beta-2) (GABPB-2) (Nuclear respiratory factor 2) (Transcription factor E4TF1-47) (Transcription factor E4TF1-53)
Other Name:	
Human Swiss Prot No.:	Q06547
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	This antibody detects endogenous levels of GABPB1 at Human, Mouse
Storage and Stability:	Store at -20° C/1 year