

Antibodies Product Datasheet

Product Name:	MAP1S Antibody
Catalog Number:	C48421
Reactivity:	Human, Mouse, Rat
Applications:	WB
Dilutions:	WB 1:500-2000
Modification:	Unmodified/Total
Source:	Polyclonal, Rabbit,IgG
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Immunogen:	Synthesized peptide derived from human MAP1S
Concentration:	1 mg/ml
Gene Name:	MAP1S BPY2IP1 C19orf5 MAP8 VCY2IP1
Protein Name:	Microtubule-associated protein 1S (MAP-1S) (BPY2-interacting protein 1) (Microtubule-associated protein 8) (Variable charge Y chromosome 2-interacting protein 1) (VCY2-interacting protein 1) (VCY2IP-1) [Cleaved into: MAP1S heavy chainMAP1S light chain]
Other Name:	
Human Swiss Prot No.:	Q66K74
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	This antibody detects endogenous levels of MAP1S at Human, Mouse,Rat
Storage and Stability:	Store at -20° C/1 year