

Antibodies Product Datasheet

Product Name:	KPNA4 Antibody
Catalog Number:	C48422
Reactivity:	Human, Mouse
Applications:	WB
Dilutions:	WB 1:500-2000
Modification:	Unmodified/Total
Source:	Polyclonal, Rabbit,IgG
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Immunogen:	Synthesized peptide derived from human KPNA4
Concentration:	1 mg/ml
Gene Name:	KPNA4 QIP1
Protein Name:	Importin subunit alpha-4 (Importin alpha Q1) (Qip1) (Karyopherin subunit alpha-4)
Other Name:	
Human Swiss Prot No.:	O00629
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	This antibody detects endogenous levels of KPNA4 at Human, Mouse
Storage and Stability:	Store at -20° C/1 year