

Antibodies Product Datasheet

Product Name:	LIPH Antibody
Catalog Number:	C48453
Reactivity:	Human, Mouse, Rat
Applications:	WB
Dilutions:	WB 1:500-2000
Modification:	Unmodified/Total
Source:	Polyclonal, Rabbit,IgG
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Immunogen:	Synthesized peptide derived from human LIPH
Concentration:	1 mg/ml
Gene Name:	LIPH LPDLR MPAPLA1 PLA1B
Protein Name:	Lipase member H (LIPH) (EC 3.1.1.-) (LPD lipase-related protein) (Membrane-associated phosphatidic acid-selective phospholipase A1-alpha) (mPA-PLA1 alpha) (Phospholipase A1 member B)
Other Name:	
Human Swiss Prot No.:	Q8WWY8
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	This antibody detects endogenous levels of LIPH at Human, Mouse,Rat
Storage and Stability:	Store at -20° C/1 year