

Antibodies Product Datasheet

Product Name:	COLEC11 Antibody
Catalog Number:	C48459
Reactivity:	Human, Mouse
Applications:	WB
Dilutions:	WB 1:500-2000
Modification:	Unmodified/Total
Source:	Polyclonal, Rabbit,IgG
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Immunogen:	Synthesized peptide derived from human COLEC11
Concentration:	1 mg/ml
Gene Name:	COLEC11 UNQ596/PRO1182
Protein Name:	Collectin-11 (Collectin kidney protein 1) (CL-K1)
Other Name:	
Human Swiss Prot No.:	Q9BWP8
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	This antibody detects endogenous levels of COLEC11 at Human, Mouse
Storage and Stability:	Store at -20° C/1 year