

Antibodies Product Datasheet

Product Name:	FZR1/Cdh1 Antibody
Catalog Number:	C48475
Reactivity:	Human, Mouse
Applications:	WB
Dilutions:	WB 1:500-2000
Modification:	Unmodified/Total
Source:	Polyclonal, Rabbit,IgG
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Immunogen:	Synthesized peptide derived from human FZR1/Cdh1
Concentration:	1 mg/ml
Gene Name:	FZR1 CDH1 FYR FZR KIAA1242
Protein Name:	Fizzy-related protein homolog (Fzr) (CDC20-like protein 1) (Cdh1/Hct1 homolog) (hCDH1)
Other Name:	
Human Swiss Prot No.:	Q9UM11
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	This antibody detects endogenous levels of FZR1/Cdh1 at Human, Mouse
Storage and Stability:	Store at -20° C/1 year