

Antibodies Product Datasheet

Product Name:	FNDC3B Antibody
Catalog Number:	C48480
Reactivity:	Human, Mouse
Applications:	WB
Dilutions:	WB 1:500-2000
Modification:	Unmodified/Total
Source:	Polyclonal, Rabbit, IgG
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Immunogen:	Synthesized peptide derived from human FNDC3B
Concentration:	1 mg/ml
Gene Name:	FNDC3B FAD104 NS5ABP37 UNQ2421/PRO4979/PRO34274
Protein Name:	Fibronectin type III domain-containing protein 3B (Factor for adipocyte differentiation 104) (HCV NS5A-binding protein 37)
Other Name:	
Human Swiss Prot No.:	Q53EP0
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	This antibody detects endogenous levels of FNDC3B at Human, Mouse
Storage and Stability:	Store at -20° C/1 year