

Antibodies Product Datasheet

Product Name:	GLRA1 Antibody
Catalog Number:	C48484
Reactivity:	Human, Rat
Applications:	WB
Dilutions:	WB 1:500-2000
Modification:	Unmodified/Total
Source:	Polyclonal, Rabbit,IgG
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Immunogen:	Synthesized peptide derived from human GLRA1
Concentration:	1 mg/ml
Gene Name:	GLRA1
Protein Name:	Glycine receptor subunit alpha-1 (Glycine receptor 48 kDa subunit) (Glycine receptor strychnine-binding subunit)
Other Name:	
Human Swiss Prot No.:	P23415
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	This antibody detects endogenous levels of GLRA1 at Human,Rat
Storage and Stability:	Store at -20° C/1 year