

Antibodies Product Datasheet

Product Name:	HLA class I (HLA-B) Antibody
Catalog Number:	C48508
Reactivity:	Human
Applications:	WB
Dilutions:	WB 1:500-2000
Modification:	Unmodified/Total
Source:	Polyclonal, Rabbit, IgG
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Immunogen:	Synthesized peptide derived from human HLA class I (HLA-B)
Concentration:	1 mg/ml
Gene Name:	HLA-B HLAB
Protein Name:	HLA class I histocompatibility antigen, B-7 alpha chain (MHC class I antigen B*7)
Other Name:	
Human Swiss Prot No.:	P01889
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	This antibody detects endogenous levels of HLA class I (HLA-B) at Human
Storage and Stability:	Store at -20° C/1 year