

Antibodies Product Datasheet

Product Name:	COG3 Antibody
Catalog Number:	C48509
Reactivity:	Human
Applications:	WB
Dilutions:	WB 1:500-2000
Modification:	Unmodified/Total
Source:	Polyclonal, Rabbit,IgG
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Immunogen:	Synthesized peptide derived from human COG3
Concentration:	1 mg/ml
Gene Name:	COG3 SEC34
Protein Name:	Conserved oligomeric Golgi complex subunit 3 (COG complex subunit 3) (Component of oligomeric Golgi complex 3) (Vesicle-docking protein SEC34 homolog) (p94)
Other Name:	
Human Swiss Prot No.:	Q96JB2
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	This antibody detects endogenous levels of COG3 at Human
Storage and Stability:	Store at -20° C/1 year