

Antibodies Product Datasheet

Product Name:	ISYNA1 Antibody
Catalog Number:	C48518
Reactivity:	Human, Mouse, Rat
Applications:	WB
Dilutions:	WB 1:500-2000
Modification:	Unmodified/Total
Source:	Polyclonal, Rabbit,IgG
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Immunogen:	Synthesized peptide derived from human ISYNA1
Concentration:	1 mg/ml
Gene Name:	ISYNA1 INO1
Protein Name:	Inositol-3-phosphate synthase 1 (IPS 1) (EC 5.5.1.4) (Myo-inositol 1-phosphate synthase) (MI-1-P synthase) (MIP synthase) (hIPS) (Myo-inositol 1-phosphate synthase A1) (hINO1)
Other Name:	
Human Swiss Prot No.:	Q9NPH2
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	This antibody detects endogenous levels of ISYNA1 at Human, Mouse,Rat
Storage and Stability:	Store at -20° C/1 year