

Antibodies Product Datasheet

Product Name:	GNAQ Antibody
Catalog Number:	C48520
Reactivity:	Human, Mouse, Rat
Applications:	WB
Dilutions:	WB 1:500-2000
Modification:	Unmodified/Total
Source:	Polyclonal, Rabbit,IgG
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Immunogen:	Synthesized peptide derived from human GNAQ
Concentration:	1 mg/ml
Gene Name:	GNAQ GAQ
Protein Name:	Guanine nucleotide-binding protein G(q) subunit alpha (Guanine nucleotide-binding protein alpha-q)
Other Name:	
Human Swiss Prot No.:	P50148
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	This antibody detects endogenous levels of GNAQ at Human, Mouse,Rat
Storage and Stability:	Store at -20° C/1 year