

Antibodies Product Datasheet

Product Name:	LASS2 Antibody
Catalog Number:	C48541
Reactivity:	Human, Mouse
Applications:	WB
Dilutions:	WB 1:500-2000
Modification:	Unmodified/Total
Source:	Polyclonal, Rabbit,IgG
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Immunogen:	Synthesized peptide derived from human LASS2
Concentration:	1 mg/ml
Gene Name:	CERS2 LASS2 TMSG1
Protein Name:	Ceramide synthase 2 (CerS2) (LAG1 longevity assurance homolog 2) (SP260) (Tumor metastasis-suppressor gene 1 protein)
Other Name:	
Human Swiss Prot No.:	Q96G23
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	This antibody detects endogenous levels of LASS2 at Human, Mouse
Storage and Stability:	Store at -20° C/1 year