

Antibodies Product Datasheet

Product Name:	Amylase Alpha Antibody
Catalog Number:	C48543
Reactivity:	Human
Applications:	WB
Dilutions:	WB 1:500-2000
Modification:	Unmodified/Total
Source:	Polyclonal, Rabbit,IgG
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Immunogen:	Synthesized peptide derived from human Amylase Alpha
Concentration:	1 mg/ml
Gene Name:	AMY2B
Protein Name:	Alpha-amylase 2B (EC 3.2.1.1) (1,4-alpha-D-glucan glucanohydrolase 2B) (Carcinoid alpha-amylase)
Other Name:	
Human Swiss Prot No.:	P19961
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	This antibody detects endogenous levels of Amylase Alpha at Human
Storage and Stability:	Store at -20° C/1 year